

T-CLASSIC

Podiaclave



The T-Classic Podiaclave has a 9 litre capacity, with a 134 °C sterilizing temperature. The internal furniture consists of a basket as standard and the unit also comes complete with a printer. The portability of T-Classic Podiaclave means that it is ideal for those practitioners who domiciliary work and need to take a sterilizer perform with them.



Model Number	T-Classic Podiaclave
Maximum Width Overall	340mm
Maximum Height Overall	335mm
Maximum net weight	6.40 kg
Chamber Diameter/Height	210mm/230mm
Chamber capacity	9 Litres
Maximum Instrument length	228mm
Maximum Load weight	3kg
Sterilisation Temperature	134 °C
Operating Pressure	2.2 Bar
Minimum Sterilisation time	3 minutes
Cycle Time Normal	20 minutes
Voltage	230 Volts
Wattage	1500 Watts
Frequency	50 – 60 Hertz
Printer	Yes
Temperature and Pressure Gauges	N/A
Furniture	Basket



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel

Model	Inner Chamber Dimensions ø x D (mm)	Volume (Litres)	Shipping Weight (kg)	Shipping Dimensions WxDxH (mm)	Single Phase Power Supply
T-Classic Podiaclave	210 x 330	9	9.25	400 x 400 x 620	220-240V, 50-60Hz 1500W, 10A

Prestige Medical Limited

East House, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

Email Sales: sales@prestigemmedical.co.uk

Email Customer Services :customerservice@prestigemmedical.co.uk

Main: +44 (0) 1254 682 622

Customer services: +44 (0) 1254 844 116

Fax: +44 (0) 1254 682 606

Sales: +44 (0) 1254 682 622